



Inaugural MyNM Conference 2018

Malaysian Society Of Nuclear Medicine & Molecular Imaging (MSNMMI)

27 – 29 September 2018 | SUNMED CONVENTION CENTRE



Singapore
General Hospital
SingHealth

Advances in detector technology in nuclear medicine & hybrid imaging

27 September 2018 Thursday (Day 1) 1430 – 1500 Hrs

S . Somanesan
Senior Principal Medical Physicist,
Dept. of Nuclear Medicine & Molecular Imaging,
Singapore General Hospital.

Hospital Radiation Safety Officer,
Chief, Radiation Response Team.
Operational & QA Manager, PET/CT Facility,



SingHealth Academic Healthcare Cluster



Singapore
General Hospital



KK Women's and
Children's Hospital



Sengkang Health



National Cancer
Centre Singapore



National Dental
Centre Singapore



National Heart
Centre Singapore



National
Neuroscience Institute



Singapore National
Eye Centre



Polyclinics
SingHealth



Bright Vision
Hospital

PATIENTS. AT THE HEART OF ALL WE DO.®



Partner in
Academic Medicine

DUKE UNIVERSITY
GRADUATE MEDICAL SCHOOL SINGAPORE

Advances of detector technology in nuclear medicine - *Scope*

- Introduction to Nuclear Medicine and PET detectors,
- Development and current trends in SPECT Detection Instrumentation:
 - Scintillators,
 - Photon Transducers
 - Semiconductors Detectors
- Summary: What the future holds for detector technology,

CZT



SiPM photo detector

Historical perspective - Developments in Instrumentation




- In the early days of nuclear medicine (1927), scans were performed by manually positioning a simple GM counter above the organ of interest.
- A major breakthrough followed with the invention of the Anger scintillation camera (early 1950's).
- The basic design of a large-area NaI(Tl) scintillation crystal, an array of photomultiplier tubes (PMTs), in combination with an absorptive collimator has been referred to as the 'Anger Camera'.

Types of gamma camera



Components of a gamma camera

- Collimator
- Sodium Iodide (NaI) crystal
- Photomultiplier tube (PMT)
- PHA circuit, discriminator & analyser circuits
- Computers for data processing & display



Development in SPECT Detection

1. Scintillators

- The advances in SPECT are directly tied to:
 - improving the quality and accuracy of the acquired projections either through:
 - better instrumentation or better correction algorithms.
- CsI(Tl), CsI(Na)[1] & LaBr₃ [2] are used in comparison with NaI(Tl).


Scintillator	Atomic number Z effective	Density ρ (g/cm ³)	Decay time (ns)	Wavelength (nm)	Relative light output (% of NaI(Tl))	Remarks
NaI (Tl)	50	3.67	200	415	100	
CsI (Tl)	54	4.5	1,000	550	45 (118*)	Long wavelength; Cardius.
CsI (Na)	54	4.51	630	420	85	Used in LinoView, s animal
LaBr ₃ :Ce	47	5.3	25	360	160	Used in s animal SPECT

- Comparable effective atomic #'s & densities \Rightarrow similar photopeak efficiency,
- Advantage: less dead edges \Rightarrow better spatial resolution.

1. Walrand S, Jamar F, de Jong M, Pauwels S. Evaluation of novel whole-body high-resolution rodent SPECT (Linoview) based on direct acquisition of linogram projections. J Nucl Med. 2005;46:1872–1880.

2. Moses WW, Shah KS. Potential for RbGd₂Br₇:Ce, LaBr₃:Ce, LaBr₃:Ce, and LuI₃:Ce in nuclear medical imaging.


Nucl Instrum Methods Phys Res A. 2005; 537:317–320.



Development in SPECT Detection

2. Photon Transducers

- PMTs have been in use for 40 years:
 - PMTs have very large electronic gains (10^6) with relatively low noise,
 - conversion efficiency is low (20%), that affects both energy resolution & intrinsic spatial resolution,
 - are affected by environmental changes such as temperature, humidity and magnetic fields and their age,
 - they are bulky as well as expensive.
- Improved replacements for conventional PMTs, [1,2]:
 - B the position sensitive PMT (PSPMT),
 - C avalanche photodiode (APD),
 - D silicon PMT [3]
- These have been recently introduced and are being incorporated into SPECT & PET detectors.





Adapted from [1]

1. Pichler BJ, Ziegler SI. Photodetectors. In: Wernick MN, Aarsvold JN, eds. Emission Tomography: The Fundamentals of SPECT and PET. 1st ed. San Diego, CA: Elsevier; 2004:255–269.
2. Bernd J. Pichler et al, Latest Advances in Molecular Imaging Instrumentation, JNM, Vol. 49 No. 6 (Suppl) June 2008
3. Moehrs S, Del Guerra A, Herbert DJ, Mandelkern MA. A detector head design for small-animal PET with silicon photomultipliers (SiPM). Phys Med Biol. 2006;51:1113–1127.

Development in SPECT Detection


3. Semiconductor Detector

- Solid-state devices (CZT) that provide direct conversion of absorbed γ -ray energy into an electronic signal [1].
- Are compact & operate at low voltages as no need for a PMT like high-gain amplification,
- The absorbed energy from a γ -ray interaction liberates charge carriers (electrons and holes) within the semiconductor.
- The induced charge on the terminals generates an electronic pulse with an amplitude proportional to the absorbed energy.
- The precision of the γ signal detection is better as no need for *light-electrical signal* conversion,
 - \Rightarrow improved energy resolution.




Advances in semiconductor detectors


- The most popular alternative for *Crystal/ PMT* is the semiconductor.
- By avoiding the need for bulky PMTs, the imaging systems can be more compact.
- The precision of the γ signal is better as no need for light-electrical signal conversion, \Rightarrow improved energy resolution.
- The most widely used semiconductors for nuclear medicine are CdTe and CdZnTe (CZT)
- The cost differential has been coming down over the decades.



Semiconductors Detectors - CdZnTe





Detector face






To cooling fin

Electronics

Gamma detection with scintillation v.s. semiconductor detectors



Scintillation Detector



Semiconductor Detector


New Approaches in SPECT

- A number of imaging systems have been introduced, that were designed specifically for cardiac imaging applications
 - With improved sensitivity in order to reduce scan times (@ 5 mins MPI).
 - Better patient comfort (imaging in an upright or reclining position)
 - By eliminating positioning the patient's arms above the head
 - Limiting claustrophobic effects (with open & smaller gantries)
 - overall sizes reduced (minimal floor space requirement).
- **Cardius XCT™**, manufactured by Digirad, Inc. Poway, CA
- **D-SPECT™**, Manufactured by Spectrum Dynamics, Haifa, Israel.




Approaches in SPECT – Alcyone™

- Manufactured by GE Healthcare using Alcyone detector
- The detectors are composed of arrays of pixilated CZT modules arranged in 180° geometry and utilising focused collimation to optimise data acquisition from the myocardium.
- Alcyone technology allows for data acquisition without detector motion.
- Spatial resolutions of 5 mm are obtainable with acquisition times of 3 to 5 minutes in typical clinical imaging situations.



PET-CT

Time Magazine's Medical Invention of the Year 2000



TIME magazine December, 2000

<http://www.prnewswire.com/news-releases/time-magazine-selects.html>



Brief Introduction of PET

- All PET systems utilize coincidence detection of the annihilation photons with 511 keV from positron decay.
- Since the paired gamma rays are anti-parallel, the detection of the gamma rays determines a line of response (LOR) along which the annihilation took place.
- The coincidence is determined by imposing an acceptance time window on the difference in time between detection of the two events.
- If the system has very high timing resolution (<1ns), ► a time-of-flight (TOF) resolution ► provide higher signal-to-noise ratio (SNR).



PET detectors

- A crystal/PMT 'block detector' uses 4 PMTs to decode an array of crystals with various combinations of reflectors and surface treatments between the crystals.
- LSO provided designers with a series of fast, bright, dense scintillators.
 - The increased light output to decode more crystals per PMT and
 - allow the use of arrays of smaller crystals, increasing resolution.
- With improved electronics and better performing PMTs, **TOF PET** appeared as a commercial product around 10 years back.



Advances in photodetectors - APD

- The “workhorse” for a scintillation detectors has been PMT. They are very high gain (typically $\sim x 10^6$), low noise, fast response and relatively moderate cost.
- The desire to further reduce the **volume** and to have detectors that can **operate within magnetic fields** has led to solid-state photodetectors.
- One of the most popular choice is the avalanche photodiode (APD).
 - It has low gain ($\sim 10^2 - 10^3$)
 - timing, suited for **non-TOF** PET applications
 - with a temperature sensitive gain (@1 or 2° celsius change)
 - SNR of current photodiode devices is a potential limitation

Advances in photodetectors - SiPM

- Silicon photomultipliers, (SiPM) are proving to be promising.
- Series of APD micro-cells.
- Gains range from 10^5 to 10^7 and have achieved 250-400 ps timing resolution with high light output, fast scintillators.
- Like most semiconductor devices, SiPMs are susceptible to **thermal noise** but noise is negligible due to high gain.
- Limitation: difficulty in scaling up to large detection areas.
- The devices can be used in magnetic fields and samples have been tested in fields up to 15 Tesla without showing any degradation.
- Can be used in TOF PET-CT and PET-MR scanners.

S. Dolinsky, et al., Nuclear Instruments & Methods in Physics Research : Timing resolution performance comparison of different SiPM devices .


Comparing PMT, APD and SiPM

- The table below compares the main features of PMT, APD and SiPM.

	PMT Photomultiplier Tube	APD Avalanche Photodiode	SiPM Silicon Photomultiplier
Quantum Efficiency	25% ... 40%	... 80%	... 80%
Compactness	Bulky	Compact	Compact
Operation Voltage	1000 – 3000 V	100 - 500 V	20 - 80 V
Gain	$10^4 - 10^9$	30 - 300	$10^5 - 10^7$
Insensitivity to Magnetic Field	No	•	•
Timing Resolution	550ps	~2000ps	<400ps
Production Costs	Medium	Low	Low

[1] D. Renker, Nuclear Instruments & Methods A 567 (2006) 48. [2] J.S. Karp, et al., Journal of Nuclear Medicine 49 (3) (2008) 462.
 [3] M. Moszynski, et al., IEEE Transactions on Nuclear Science NS53 (5) (2006) 2484; [4] A Iltis, et al., Nuclear Instruments & Methods
 A 563 (2006) 359. [5] C.L. Melcher, Journal of Nuclear Medicine 41 (6) (2000) 1051.

Advantages of the Silicon Photomultiplier



- Single photon resolution
- High photon detection efficiency
- High gain and high signal to noise ratio
- Very fast device with short recovery and good timing
- Low bias voltage
- Insensitive to magnetic fields
- Very low nuclear counting effects
- Very compact devices
- Simple calibration and monitoring
- Not damaged by day light
- Low costs



[1] D. Renker, Nuclear Instruments & Methods A 567 (2006) 48. [2] J.S. Karp, et al., Journal of Nuclear Medicine 49 (3) (2008) 462. [3] M. Moszynski, et al., IEEE Transactions on Nuclear Science NS53 (5) (2006) 2484; [4] A Iltis, et al., Nuclear Instruments & Methods A 563 (2006) 359. [5] C.L. Melcher, Journal of Nuclear Medicine 41 (6) (2000) 1051.

So What in the Market now for hybrid PET-MR & PET-CT scanners?



Philips™ Vereos™ PET/CT


- Using Digital Photon Counting (DPC) technology so every scintillator crystal is directly coupled to a single SiPM detector.
- Each detector light sensor contains microcells that are capable of detecting a single photon converting the photon signal directly to an electrical output.
- The system contain more than 400,000 cells.
- Fast timing resolution of 345ps so the TOF localization accuracy is 5.2cm.



Courtesy of Philips

Siemens™ Biograph™ mMR

- First fully integrated PET/MR in the world.
- Currently uses LSO crystal with APDs for the PET/MR system.
- Time of flight is not available due to the limitation of APDs.



Images Courtesy of Siemens Medical Systems

GE™ SIGNA™ PET/MR

- Simultaneous time of flight (TOF) PETMR
- Currently uses LySO crystal with Silicon Photomultiplier (SiPM)



PET/MR Module Design

Light Tight RF Shield with Copper Coating



25mm Scintillator (LBS) Crystal Array with Light Guides and Enhanced Spectral Reflectors (ESR)



SiPM – Silicon Photomultiplier with Circuit Boards/ASICs



Thermal Coupling Gasket



Aluminum Mounting Plate



MR compatible casing with copper coating



Ceramic Cold Plate with copper coating



Distributed Module Electronics



Images Courtesy of GE Healthcare




Inaugural MyNM Conference 2018

Malaysian Society Of Nuclear Medicine & Molecular Imaging (MSNMM)

27 – 29 September 2018 | SUNMED CONVENTION CENTRE


Summary: Future of photodetectors

- SiPMs have shown good energy and timing resolution, and an ability to decode arrays of scintillator elements,
- SiPM's have confirmed their potential as possible replacements for PMTs in PET detectors.
- Key questions that will determine whether PET scanners of the future will continue to be based on SiPMs:
 - Can SiPMs arrays be made with larger sensitive areas of at least $20 \times 20 \text{ mm}^2$ without any degradation;
 - Can performance and cost per unit area be significantly less than PMT's?
 - If the answers are 'Yes', then it would be hard to imagine that we will not see widespread adoption of SiPM technology in PET in the future



Summary: Future of solid state detectors

- This technology was approved by the FDA since 1997.
- The high cost and sporadic availability of CZT crystals have limited the adoption by the majority of companies & use in large FOV scanners.
 - Currently a vendor has a whole body dual detector CZT scanner for general SPECT and whole body imaging.
- One strong driver for solid-state technology could be its ability to be coupled with other imaging modalities, such as MRI,
- In addition, solid-state detectors have good low-energy resolution.
- **While the ideal detector has yet to be realized, the capabilities of modern scanners have continued to advance.**



Reference:

- Tom K Lewellen, *Recent developments in PET detector technology*, *Phys Med Biol.* 2008 Sep 7; 53(17): R287-R317.
- Ron Jaszczak et al. *The early years of single photon emission computed tomography: an anthology of selected reminiscences*. *Phys.Med.Biol.* 51(2006) R99-R115
- Tom K Lewellen. *Recent developments in PET detector technology*. *Phys Med Biol.* 2008 september 7;53(17):R287-R317
- J Huizenga et al. *A fast preamplifier concept for SiPM-based time-of-flight PET detectors*. *Nuclear Instruments and Methods in Physics Research A* 695 (2012)379-384
- Mark T.Madsen. *Recent advances in SPECT Imaging*. *J Nucl Med* 2007,48:661-673
- Todd E Peterson et al. *SPECT detectors: the Anger Camera and beyond*. *Phys. Med. Biol* 56(2011) R145-R182
- Emilie Roncali & Simon R. Cherry; *Annals of Biomedical Engineering*, Vol. 39, No. 4, April 2011 (2011) pp. 1358–1377
- D. Renker, *Nuclear Instruments & Methods A* 567 (2006) 48.
- J.S. Karp, et al., *Journal of Nuclear Medicine* 49 (3) (2008) 462.
- M. Moszynski, et al., *IEEE Transactions on Nuclear Science* NS53 (5) (2006) 2484;
- A Iltis, et al., *Nuclear Instruments & Methods A* 563 (2006) 359.
- C.L. Melcher, *Journal of Nuclear Medicine* 41 (6) (2000) 1051.
- Walrand S, Jamar F, de Jong M, Pauwels S. *Evaluation of novel whole-body high-resolution rodent SPECT (Linoview) based on direct acquisition of linogram projections*. *J Nucl Med.* 2005;46:1872–1880.
- Moses WW, Shah KS. *Potential for RbGd₂Br₇:Ce, LaBr₃:Ce, LaBr₃:Ce, and LuI₃:Ce in nuclear medical imaging*.
- *Nucl Instrum Methods Phys Res A.* 2005; 537:317–320.





Inaugural MyNM Conference 2018

Malaysian Society Of Nuclear Medicine & Molecular Imaging (MSNMMI)

27 – 29 September 2018 | SUNMED CONVENTION CENTRE

