

All About Hanau H2 And WideVue Articulators – PCC 25

Chief Mentor

Dr. Anjali Bhoyar

Co-Mentor

Dr.SeemaSathe

Suitable for- New Learners & Intermediate learners

Learning objectives

The participants shall be able to learn -

1. Functions, requirements, limitations and types of articulators
2. Parts of Semi - adjustable Articulators Viz. Hanau H2 and WideVue Model
3. Articulator Programming and Operation
4. Cleaning, Disinfection and Maintenance of Articulators

Background information – Articulators are the mechanical devices used in dentistry, especially in Prosthetic dentistry for Diagnosis, treatment planning and fabrication of artificial restorations or prosthesis. They act as patients in absence of patients, thus saving patient's time and operator's chairside time. The beginners shall be well acquainted with the device in order to operate them correctly.

Important topics/concepts covered - This preconference Course throws light on Two

of the commonly used articulators in prosthetic dentistry such as HANAU H2 AND WIDEVUE SEMIADJUSTABLE ARTICULATORS.

Whether only lecture or includes any video demonstration – Includes Video Demonstration

Any instructions for the participants

Participants are advised to carry Hanau H2 / WideVue Articulators with them