


Why Mint Lift® is so effective?


Mint Lift® is made of **safe and absorbable PDO material**.

During the treatment, **Mint Lift® 360°, 3D helical structure cogs threads** will be implanted into the **SMAS layer** (Superficial Musculo-Aponeurotic System) to **instant lifting** your face upwards; and to **stimulate collagen** to make the effect **long-lasting**.



Who is the target of Mint Lift® ?

- Those who want to improve wrinkles and sagging look
- Want a V face line and a tighter jawline
- Those who are worried about the risk of surgery
- Those who want immediate lifting effect



How long will it last? How long does it take to absorb?

After Mint Lift® is implanted, it will be **gradually absorbed by the human body** in about 8-10 months, and the effect can last for about 1.5-2 years.

However, it **depends** on individual's situation.

It's recommended to do a follow up after six months.



Will it hurt? How long is the treatment downtime?

Mint Lift® is a **new generation of minimally invasive lifting**. The wound is very small. The range of pain **can be tolerated** by all.

Minimal to no downtime, so you can wear makeup right after the treatment.



MINT®

Minimally Invasive Non-surgical Thread

1st Thread

Received **3 FDA Cleared**
For Lifting & Wrinkles Improving



Brows • Forehead • Double Chin
Nasolabial Folds • Dropping cheek meat



U.S. FDA
3 Cleared



Natural &
Long Lasting Result



Lifting & Stimulating
Collagen



Non-Surgical
Treatment



30 mins
Treatment/Part




CE Approval

HANS BIOMED HansBiomed Corp
Leading Regenerative Medicine

Orient EuroPharma (M)

myoeaesthetic

MINT™ EASY



TAPER
POINT



Strong
Traction

Combined
with
Needle

Strength	Technical Level	Design Diversity
----------	-----------------	------------------



Thread Length	USP Size	Needle Length	Needle Gauge	Package
---------------	----------	---------------	--------------	---------


424mm

1

130mm

19G(Needle)

2ea*1Pack/Box



Cog Design: Bi-direction

01 No.1 Best Brand

Customer satisfaction awarded
for 3 consecutive years!



02 Lifting & Wrinkles Improving Approval

The approval including
U.S. FDA 510(k) 3 Cleared and CE 2265.



UNITED STATES
K130191, K192423 &
K220549
3 Cleared


CE
2265

MINT® Minimally Invasive Non-surgical Thread



03 Patented

Patent manufacturing process
No.: CPCTP11-0007.



Press Molding Method

04 Outstanding Strength

4 times stronger than other
conventional cutting threads.




Other Brands



MINT Lift®

MINT™ FINE



W Type



Delicate
Areas
Lifting

Strength	Technical Level	Design Diversity
----------	-----------------	------------------



Thread Length	USP Size	Cannula Length	Cannula Gauge	Package
---------------	----------	----------------	---------------	---------


150mm

2-0

100mm

18G

4ea*5Packs/Box



Cog Design: Bi-direction

05 360°, 3D Helical Structure Cogs

The 3D barbs are molded and positioned
helically in 360° on the thread surface
to achieve a long-lasting and powerful
lifting effect.



Other Brands

MINT Lift®



360°

